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

013429RE-1CP Reference electrode (Ag/AgCl/Saturated KCl)


1 Contents of electrode

- a) RE-1CP main body
- b) Holder cover
- c) Protective cap
- d) Holder  
(porous glass liquid junction)

2. How to use


2.1 Remove the protective cap and check air bubble in the holder. If air bubbles were present around the porous glass, slightly flick the electrode to clear the bubbles.

2.2 If the inner solution doesn't reach the main body, inject the KCl saturated solution until the middle position of upper of holder.

2.3 Expose the side hole of holder (d) in measuring.

3 How to keep

After use, wash the tip of electrode by ion exchange water. Cover the side hole of holder with holder cover and soak in KCl saturated solution. If don't use for the long time, disassemble the electrode and sink the main body in KCl saturated solution and sink the holder filled with ion exchange water in ion exchange water immerse the tip of holder in ion exchange water.

012108 RE-PV Preservative vial for reference electrode (option)

!!CAUTION!!

1. Do not remove the label of electrode. The label is required for any case of inquiry.
2. Refrain the electrode from using in strong acid or alkaline solution.
3. Use the electrode at room temperature and atmospheric pressure.
4. To avoid breaking the glass tube, the electrode should be protected from strong shock.
5. The porous glass of tip may be discolored after use, which is due to the nature of the glass. Hence, we can not respond to exchange.
6. This electrode is intended for aqueous. Do not use it in organic solvent.

You can browse the checking data of electrode in below URL.  
<http://www.als-japan.com/dl/index.htm>