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

013458 RE-2BP Calomel Reference electrode

1 .Contents of electrode


- a) RE-2BP main body
- b) Holder cover
- c) RE-2BP holder
- d) Protective cap

2. Check the electrode before use


2.1 Check whether there are air bubbles in the solution of RE-2BP holder(c) firstly. If there are air bubbles around the porous glass or the tip of the RE-2BP main body(a), slightly flick the electrode to remove the bubbles.

2.2 Please expose the side hole of holder(c) by lowering the holder cover during the measurement.

Note : In a case of the RE-2BP main body(a) being exposed due to deficit of internal solution after using, let the holder cover (b) lower and inject a saturated KCl solution through the side hole until the RE-2BP main body(a) is soaked in solution.

3 .How to keep

After use, wash the tip of the RE-2BP holder(c) with ion-exchanged water. Cover the side hole of holder with holder cover(b) and soak the electrode in saturated KCl solution.

If don't use the electrode for long time, please seal the connecting part between RE-2BP main body(a) and RE-2BP holder(c) with parafilm and soak in saturated KCl solution. Please take care of possibility that inadequate storing may cause a shift of electrode potential and damage of liquid junction.

Keep in 012108
Preservative vial (option)

!!CAUTION!!

1. The electrode consists of hazardous substance. Be careful handling.
2. Do not remove the label of electrode. The label is required for any case of inquiry.
3. Refrain the electrode from using in strong acid or alkaline solution.
4. Use the electrode at room and under atmospheric temperature and atmospheric pressure.
5. To avoid breaking, the electrode should be protected from strong shock.
6. Don't disassemble the electrode. If disassemble the electrode, we cannot warranty.
7. The holder tip may be discolored after use, which is due to the nature of the ceramics. Hence, we can not respond to exchange.
8. This electrode is intended for aqueous. Do not use it in organic solvent.

You can browse the checking data of electrode in below URL.
<http://www.als-japan.com/d/>

Manufactured by: ALS Co., Ltd Exclusive distributor: BAS Inc.
URL: <http://www.als-japan.com>
E-mail: sales@als-japan.com