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


1. The structure of electrode

- a. RE-61AP main body
 b. Holder cover
 c. Protective cap
 d. Holder
 (ceramics liquid-junction)
 The inner solution is the ion-exchange water at time of receipt. **Please prepare inner solution (1 M NaOH) before use.**

2. How to use


2.1 When you receive the RE-61AP, disassemble the electrode and discard water in the holder. Then wash the Holder (d) with a few portions of inner solution (1 M NaOH) three times, then inject the inner solution to be 5mm below the level of side hole. If you see air bubbles in the holder, remove the bubbles by flicking or shaking.


2.2 Fix the main body (a) to the holder (d) tightly. Remove a little of internal solution if it seems overflow from side hole of the holder. (Note: the side hole is for electrode inner pressure adjustment.)

2.3 Cover the side hole of holder (d) with the holder cover (b), soak the electrode in the solution identical with inner solution to stabilize the potential. It will take several days until potential is stabilized.

2.4 The stability of the electrode potential can be checked using potentiometry against the standard reference electrode (SCE etc.) in same solution of inner solution. The indication of the potentiometry against the SCE is $-115 \text{ mV} \pm 20 \text{ mV}$ (25 °C) with using 1 M NaOH.

2.5 Please use the electrode after checking the stability of potential. Expose the side hole of holder (d) in measuring.

3. How to keep

After use, wash the tip of electrode with ion-exchange water. Cover the side hole of holder with holder cover and soak in 1 M NaOH solution. When you don't use it for long period, please disassemble the electrode, keep the main body in 1M NaOH, fill ion-exchange water in the holder then immerse the tip of holder in an ion-exchange water.

!!CAUTION!!

1. The electrode consists of harmful agents. Be careful handling. At the time to discard the inner solution after final usage, please comply with the local government's law.
2. Do not remove the label of electrode. The production serial code is required for any case of inquiry.
3. The electrode is for alkaline solution measurement. Can't use in organic solvent.
4. Use the electrode at ordinary temperatures and pressures.
5. Avoid a strong shock to the electrode.
6. The electrode has individual difference and sensitive to the temperature, in some case the electrode potential may not coincide with theoretical value.